



Product description

Order number	95.64.924
Product description	Storage box, slip-on lid, material: cardboard, white, format: 12 x 12, for 144 collection tubes, suitable for tubes up to max. 9.5 mm Ø, tube height between 19-26 mm, 1 piece(s)/bag

Product characteristics

Closure type	slip-on lid
Suitable for	tubes up to max. 9.5 mm Ø, tube height between 19-26 mm
Storage capacity	144
Format	12 x 12

Size

Width of product	135 mm
Height of product	28 mm
Length of product	135 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	216
Piece(s) / pallet	2592
Depth of case	521 mm
Width of case	521 mm
Height of case	610 mm
Case volume	0.1653 cbm
Weight of product	0.0632 kg
Weight of case	15.6 kg
EAN of inner box	4038917095723
EAN case	4038917067461

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

