



Product description

Order number	95.64.998
Product description	Storage box, slip-on lid, material: cardboard, white

Product characteristics

Label/ Print	with print
Closure type	slip-on lid

Size

Width of product	135 mm
Height of product	45 mm
Length of product	135 mm

Material & colours

Product material	cardboard
Colour of product	white
Closure material	cardboard
Colour of cap	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Storage box, slip-on lid, cardboard

Page 2

Purity & certification

Product category	no medical device IVD
Batched	yes

Packaging

Minimum order qty.	10
Piece(s) / inner box	10
Piece (s) / outer case	120
Piece(s) / pallet	2160
Case volume	0.1134 cbm
Weight of product	0.065 kg
Weight of case	7.86 kg
EAN of inner box	4038917246927
EAN case	4038917246910

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-25

This document was prepared by EDP support and is valid without signature.



SARSTEDT