




**Product Specification**  
**Double block rack D17, Ø opening: 17 mm, 10 x 10, grey**

**Page 1**



**Product description**

|                     |  |
|---------------------|--|
| Order number        | 93.852.200   |
| Product description | Double block rack D17, for 100 collection tubes, Ø opening: 17 mm, format: 10 x 10, grey, material: PP |

**Product characteristics**

|                        |         |
|------------------------|---------|
| Number of wells        | 100     |
| Storage capacity       | 100     |
| Format                 | 10 x 10 |
| Autoclavable (max. °C) | 121 °C  |

**Size**

|                     |        |
|---------------------|--------|
| Diameter of opening | 17 mm  |
| Width of product    | 218 mm |
| Height of product   | 45 mm  |
| Length of product   | 209 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | grey               |

Purity & certification

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | no                      |

Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1             |
| Type of smallest subpackaging | -             |
| Piece(s) / inner box          | 1             |
| Piece (s) / outer case        | 12            |
| Piece(s) / pallet             | 528           |
| Depth of case                 | 453 mm        |
| Width of case                 | 433 mm        |
| Height of case                | 165 mm        |
| Case volume                   | 0.0324 cbm    |
| Weight of product             | 0.3917 kg     |
| Weight of case                | 5.41 kg       |
| EAN of inner box              | 4038917122146 |
| EAN case                      | 4038917243124 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

