



**Product Specification**  
**Protective cover, ABS, light grey, suitable for**  
**Block rack D17**

**Page 1**



**Product description**

Order number	93.857
Product description	Protective cover, material: ABS, light grey, suitable for Block rack D17, 25 piece(s)/case

**Product characteristics**

Suitable for	Block rack D17
Autoclavable (max. °C)	121 °C

**Size**

Width of product	108 mm
Length of product	208 mm

**Material & colours**

Product material	Acrylnitrile butadiene styrene (ABS)
Colour of product	light grey

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Protective cover, ABS, light grey, suitable for**  
**Block rack D17**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	25
Piece(s) / pallet	800
Depth of case	548 mm
Width of case	368 mm
Height of case	225 mm
Case volume	0.0454 cbm
Weight of product	0.0656 kg
Weight of case	2.2 kg
EAN of inner box	4038917084420
EAN case	4038917078054

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-05-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**