




**Product Specification**  
**Rack, styrofoam, format: 20 x 5, suitable for**  
**S-Monovette® Ø 11 mm**

**Page 1**



**Product description**

Order number	95.1439
Product description	Rack, material: styrofoam, white, format: 20 x 5, (LxWxH): 338 x 90 x 50 mm, suitable for S-Monovette® Ø 11 mm, 1 piece(s)/blister

**Product characteristics**

Suitable for	S-Monovette® Ø 11 mm
Format	20 x 5

**Size**

Width of product	90 mm
Height of product	50 mm
Length of product	338 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Rack, styrofoam, format: 20 x 5, suitable for**  
**S-Monovette® Ø 11 mm**

**Page 2**

**Material & colours**

Product material	styrofoam
Colour of product	white

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	40
Piece(s) / pallet	640
Depth of case	498 mm
Width of case	428 mm
Height of case	380 mm
Case volume	0.081 cbm
Weight of product	0.0235 kg
Weight of case	1.7 kg
EAN of inner box	4038917086646
EAN case	4038917005173

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**