



Product description

Order number	95.64.249
Product description	Rack, lid enclosed, material: styrofoam, white, format: 20 x 5, (LxWxH): 340 x 100 x 65 mm, for 100 collection tubes, suitable for tubes Ø 10.8 mm, 1 piece(s)/blister

Product characteristics

Suitable for	tubes Ø 10.8 mm
Storage capacity	100
Format	20 x 5

Size

Width of product	100 mm
Height of product	65 mm
Length of product	340 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	styrofoam
Colour of product	white
Closure material	styrofoam
Colour of cap	white

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	30
Piece(s) / pallet	480
Depth of case	498 mm
Width of case	428 mm
Height of case	380 mm
Case volume	0.081 cbm
Weight of product	0.0507 kg
Weight of case	2.3 kg
EAN of inner box	4038917095570
EAN case	4038917007757

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

