



Product description

Order number	95.64.257
Product description	Rack, material: styrofoam, white, format: 5 x 5, (LxWxH): 204 x 167 x 32 mm, for 25 collection tubes

Product characteristics

Storage capacity	25
Format	5 x 5

Size

Width of product	167 mm
Height of product	32 mm
Length of product	204 mm

Material & colours

Product material	styrofoam
Colour of product	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Rack, styrofoam, format: 5 x 5

Page 2

Packaging

Minimum order qty.	24
Piece(s) / inner box	24
Piece (s) / outer case	24
Piece(s) / pallet	864
Depth of case	413 mm
Width of case	464 mm
Height of case	258 mm
Case volume	0.0492 cbm
Weight of case	0.5 kg
EAN of inner box	4038917095624
EAN case	4038917007801

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT