



## Product Specification

Forensic swab, pointed, in the tube with ventilation membrane, ISO 18385, 95 mm, cotton

Page 1



### Product description

Order number	80.629.101
Product description	Forensic swab, pointed, in the tube with ventilation membrane, with paper label, ETO treated, ISO 18385, passive drying, stem length: 95 mm, stem material: wood, swab tip pointed, swab material: cotton, tube: PP, blister packaging aluminium free, 1 piece(s)/blister

### Product characteristics

Label/ Print	with paper label
Cap	Screw cap with cotton bud holder
Closure type	screw cap
Design	in the tube with ventilation membrane
Application	Removal of forensic samples

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-28

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Size

Diameter	18 mm
Width of product	18.05 mm
Height of product	106.6 mm
Length of product	107 mm
Length including cap	97.6 µm
Length excluding cap	105.4 mm
Length including cap and plunger	106.6 µm
Length of stem	95 mm

## Material & colours

Colour of product	multi-coloured
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Swab material	cotton
Stem material	wood
Tube material	Polypropylene (PP)

## Purity & certification

Satisfies the requirement	ISO 18385
Product category	no medical device   IVD
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Minimum order qty.	800
Type of smallest subpackaging	blister
Piece(s) / inner box	100
Piece (s) / outer case	800
Piece(s) / pallet	20000
Depth of box	211 mm
Width of box	151 mm
Height of box	171 mm
Depth of case	437 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Forensic swab, pointed, in the tube with  
ventilation membrane, ISO 18385, 95 mm, cotton**

**Page 3**

Width of case	363 mm
Height of case	335 mm
Case volume	0.0531 cbm
Weight of product	0.0065 kg
Weight of case	7.29 kg
EAN of inner box	4038917315579
EAN case	4038917428972

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2024-03-28

This document was prepared by EDP  
support and is valid without signature.



**SARSTEDT**