



Product description

Order number	90.2620.302
Product description	length: 1.5 m, transparent, 1 piece(s)/bag

Product characteristics

Suitable for	IVARO
--------------	-------

Size

Diameter	0.8 mm
Length of product	1.5 m

Material & colours

Product material	Bioprene
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	10
Case volume	0.0062 cbm
Weight of product	0.0005 kg
EAN of inner box	4038917584197
EAN case	4038917584180

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

