



Product Specification

Container, Sodium thiosulphate, 250 ml, (LxØ): 78 x 70 mm, graduated, PS, transparent, with paper label

Page 1



Product description

Order number	51.9923.821
Product description	Container, Sodium thiosulphate, max. working volume: 250 ml, (LxØ): 78 x 70 mm, transparent, cap: yellow, graduated, material: PS, screw cap, cap assembled, with paper label, label/print: black, sterile, 210 piece(s)/bag

Product characteristics

Type of preparation	Sodium thiosulphate solution
Label/ Print	with paper label
Colour of print/label	black
Graduation	yes
Cap	standard screw cap
Closure type	screw cap
Application	Container for water samples
Base shape	flat base
Volume of work	250 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

Container, Sodium thiosulphate, 250 ml, (LxØ): 78 x 70 mm, graduated, PS, transparent, with paper label

Page 2

Size

Sample volume	250 ml
Diameter	70 mm
Height of product	78 mm
Length of product	78 mm

Material & colours

Product material	Polystyrene (PS)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	yellow
Tube material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	sealed at 9"hg/10min
Latex-free	yes
Purity standard	sterile
Pyrogen-free	no
Safety data sheet	no

Packaging

Minimum order qty.	210
Type of smallest subpackaging	bag
Piece(s) / inner box	240
Piece (s) / outer case	210
Piece(s) / pallet	5040
Depth of case	534 mm
Width of case	384 mm
Height of case	519 mm
Case volume	0.1025 cbm
Weight of case	10.15 kg
EAN of inner box	4038917154567
EAN case	4038917235068

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-16

This document was prepared by EDP support and is valid without signature.



SARSTEDT