To make your workflow excellent.

SARSTEDT

Certificate
Cryo Performance Tested quality



For the use of SARSTEDT consumables in cryopreservation

Cryo Performance Tested items fulfil the following certified purity criteria:

- ✓ Sterile
- ✓ Pyrogen-free/endotoxin-free
- ✓ DNA-free
- ✓ DNase/RNase-free
- ✓ Non-cytotoxic
- ✓ Non-mutagenic
- ✓ CE IVD

We guarantee that the following limit values are met:

Sterility validated according to ISO 11137
Pyrogens/endotoxins <0.06 EU/ml
Human DNA <0.5 pg/µl, bacterial DNA <0.02 pg/µl
DNase <$1x10^{-5}$ U/µl, RNase <$1x10^{-9}$ Kunitz-units/µl
Non-cytotoxic according to ISO 10993
Non-mutagenic according to the Ames Test II

SARSTEDT AG & Co. KG

Central Quality Assurance

Technical modifications reserved

GDO_20_008_0100_NA_NF